


*HA coated screw-cup system for conical taper-fitting ceramic or PE liners.*



**Material:**

**Acetabular cup :** HA coated alumina grit-blasted titanium TA6V alloy.

**Liners :**

- UHMW Polyethylene
- Alumina ceramic.

**Packaging :** vacuum packed and gamma ray sterilized.

Manufactured by Evolutis S.A.S  
Avenue de la Libération - F- 42720 BRIENNON

Contact SCANOS at

☎ : +33.5.61.73.67.16 📠 : +33.5.61.73.52.35

e-mail : [biomed.biomaterials@wanadoo.fr](mailto:biomed.biomaterials@wanadoo.fr)

visit us at : [www.scanos.com](http://www.scanos.com)



**Distributed by:**

### CAPTIV SCW HA SCREW-CUP

Order ref	Description
H29 0044	Cup diameter 44 mm
H29 0046	Cup diameter 46 mm
H29 0048	Cup diameter 48 mm
H29 0050	Cup diameter 50 mm
H29 0052	Cup diameter 52 mm
H29 0054	Cup diameter 54 mm
H29 0056	Cup diameter 56 mm
H29 0058	Cup diameter 58 mm
H29 0060	Cup diameter 60 mm
H29 0063	Cup diameter 63 mm
H29 0066	Cup diameter 66 mm
H29 0069	Cup diameter 69 mm


**A rational easy to use clear layout out instrument set to gain time both operatively and with cleansing.**

### POLYETHYLENE LINERS

Order ref	Description
H29 2244	Spherical liner for a 22.2 head and a 44 mm cup
H29 2246	Liner for a 22.2 head and a 46 mm cup
H29 2848	Liner for a 28 head and a 48 mm cup
H29 2850	Liner for a 28 head and a 50 mm cup
H29 2852	Liner for a 28 head and a 52 or 54 mm cup
H29 2856	Liner for a 28 head and a 56 or 58 mm cup
H29 2860	Liner for a 28 head and a 60 or 63 mm cup
H29 2866	Liner for a 28 head and a 66 or 69 mm cup

### ALUMINA CERAMIC INLAYS


Order ref	Product description
H29 A2848	Alumina ceramic inlay for 28 mm head for cup 48 mm
H29 A2850	No longer available. Only H29 A3250 now used.
H29 A2852	No longer available. Only H29 A3252 now used
H29 A3250	Alumina ceramic inlay for 32 mm head for cup 50
H29 A3252	Alumina ceramic inlay for 32 mm head for cup 52 or 54 mm
H29 A3256	Alumina ceramic inlay for 32 mm head for cup 56 or 58 mm
H29 A3260	Alumina ceramic inlay for 32 mm head for cup 60 or 63 mm
H29 A3266	Alumina ceramic inlay for 32 mm head for cup 66 or 69 mm



### INSTRUMENT SET

Order ref	Description
H29 9000	CAPTIV autoclave case
S01 002	CAPTIV autoclave case lid
H30 002	Alumina liner insertion device
H30 001	Alumina liner extractor
H44 001	Cup ratchet handle
H03 0300	Click-on adaptor
H03 0328	liner impactor head
H03 0422	Ø22.2 centraliser
H03 0428	Ø28 centraliser
H03 0432	Ø32 centraliser
S01 005	3.5 mm hex screwdriver

**IMPORTANT:** Prior to inserting the Alumina inlay check that the inside of the cup is clean. The use of the suction cup instrument is essential in order to ensure a correct fit. To be able to remove the device without compromising the fixation of the inlay pull trigger to release suction vacuum.



Manufactured by Evolutis S.A.S  
Avenue de la Libération - F- 42720 BRIENNON



Contact SCANOS at  
☎ : +33.5.61.73.67.16 ☎ : +33.5.61.73.52.35  
e-mail : [biomed.biomaterials@wanadoo.fr](mailto:biomed.biomaterials@wanadoo.fr)  
visit us at : [www.scanos.com](http://www.scanos.com)

A product from...

**SCANOS**