

A product from...

**SCANOS**

# CAPTIV' ELS

## Double split cup

The CAPTIV' ELS cup is an 'elastic' double split cup that ensures :

- an **exact apical dome fit** : where anatomical constraints are most concentrated.
- an upper and peripheral **press-fit** : where primary stability is most essential.
- elasticity **between the two horns**: where the acetabulum is naturally flexible.

The CAPTIV' ELS cup is a non-cemented exact-fit press-fit cup. High secondary stability is ensured by a strong 100µm uniform HA full-coating on a 20µm rugosity porous titanium metal back substrate surface.



Poor bone quality in some patients may require the use of bone screws in order to enhance primary stability and ensure placement.

### Material :

- Cup: TA6V ELI titanium alloy in compliance to ISO 5832-3 standard.
- Inlay: UHMW polyethylene in compliance to ISO 5834-1 & 2 standards.
- Screws: TA6V ELI titanium alloy in compliance to ISO 5832-3 standard.

**Packaging** : vacuum packed and gamma ray sterilized.

### Distributed by:



Manufactured by:


**evolutis**  
42720 BRIENNON

## Double split cup

Order ref.	Description
H59 0044	Cup size 44 mm
H59 0046	Cup size 46 mm
H59 0048	Cup size 48 mm
H59 0050	Cup size 50 mm
H59 0052	Cup size 52 mm
H59 0054	Cup size 54 mm
H59 0056	Cup size 56 mm
H59 0058	Cup size 58 mm
H59 0060	Cup size 60 mm
H59 0062	Cup size 62 mm
H59 0064	Cup size 64 mm

## POLYETHYLENE inlays: flat or peaked.

Description	Spherical Ø 22.2	Flat Ø 28	Peaked Ø 28
Inlay for cup size Ø44 mm	H59 P2244		
Inlay for cup size Ø46 mm	H59 P2246		
Inlay for cup size Ø48 mm	H59 P2248	H59 P2848	H59 2848
Inlay for cup size Ø50 mm	H59 P2250	H59 P2850	H59 2850
Inlay for cup size Ø52 mm	H59 P2252	H59 P2852	H59 2852
Inlay for cup size Ø54 mm	H59 P2254	H59 P2854	H59 2854
Inlay for cup size Ø56 mm	H59 P2256	H59 P2856	H59 2856
Inlay for cup size Ø58 mm	H59 P2258	H59 P2858	H59 2858
Inlay for cup size Ø60 mm	H59 P2260	H59 P2860	H59 2860
Inlay for cup size Ø62 mm	H59 P2262	H59 P2862	H59 2862
Inlay for cup size Ø64 mm	H59 P2264	H59 P2864	H59 2864



## CUP SCREWS (sterile, packed individually)

Order ref.	Description
H15 S6020	Sterile screw Ø 6 length 20 mm
H15 S6025	Sterile screw Ø 6 length 25 mm
H15 S6030	Sterile screw Ø 6 length 30 mm
H15 S6035	Sterile screw Ø 6 length 35 mm
H15 S6040	Sterile screw Ø 6 length 40 mm
H15 S6045	Sterile screw Ø 6 length 45 mm

## INSTRUMENT SET

Order ref:	Description
H30 9001	CAPTIV PRF-ELS autoclave
H03 02XX	Metal frame shell gauge 44-60mm
H03 026	Cup holder and impactor
H03 007	Angle indicator
H03 015	Peaked liner impaction plate
H03 016	Flat liner impaction plate
H03 023	Ø 22.2 impaction head
H03 024	Ø 28 impaction head
H03 025	Ø 32 impaction head
S01 001	Flexible AO adaptor
S01 002	Ø 3,2 mm drill: length 35 mm
S01 010	Ø 3,2 mm drill: length 25 mm
S01 004	Ø 6 mm drill: length 25 mm (not used)
H03 004	Screw hole depth-gauge.
H34 5732	Screw and stud drill guide
S01 007	Articulated 3.5 mm hex screwdriver
S01 005	M8 AO impactor adaptor
H30 001	Bilox inlay extractor (not used)
H30 016	Bilox inlay suction inserter (not used)