


# CAPTIV PRF

A product from...

## SCANOS

*Exact-fit press-fit peripheral cup with a conical direct fit for alumina ceramic inlays.*



A multiple outer shell and liner option cup

Fluide tight screw holes and apical dome hole with porous sterile titanium plugs.



150µm hydroxyapatite coating on a 20µm porous titanium sub-layer surface.



Threaded dome screw holes available for proper screw placement and apical dome hole option to engage proper seating depth.

**For a perfect press-fit, use the same size cup as the last size reamer used.**

**Material:**

**Cup, screws and plugs:** titanium TA6V ELI alloy in compliance ISO 5832-3 standards,  
Cup coating: HAC coated.

**Inlays :** Alumina ceramic in compliance to NFS 90-408 standards or UHMW PE in compliance to ISO 5834 -1 & 2 standards

**Packaging:** vacuum packed and gamma ray sterilized.

A SCANOS product

☎: +33.5.61.73.67.16 ☎ : +33.5.61.73.52.35

e-mail : [biomed.biomaterials@wanadoo.fr](mailto:biomed.biomaterials@wanadoo.fr)

visit us at : [www.scanos.com](http://www.scanos.com)

**CE**  
**0499**

Manufactured by Evolutis S.A.S  
Avenue de la Libération - F- 42720 BRIENNON

## HAC COATED CUP

Order ref.	Description
H29 PF44	Cup size 44 mm
H29 PF46	Cup size 46 mm
H29 PF48	Cup size 48 mm
H29 PF50	Cup size 50 mm
H29 PF52	Cup size 52 mm
H29 PF54	Cup size 54 mm
H29 PF56	Cup size 56 mm
H29 PF58	Cup size 58 mm
H29 PF60	Cup size 60 mm
H29 PF63	Cup size 63 mm
H29 PF66	Cup size 66 mm
H29 PF69	Cup size 69 mm

## SCREWS - PEGS - PLUGS

Order ref.	Description
H15 SA6020	Ø 6 mm screw :20 mm
H15 SA6025	Ø 6 mm screw :25 mm
H15 SA6030	Ø 6 mm screw :30 mm
H15 SA6035	Ø 6 mm screw :35 mm
H15 SA6040	Ø 6 mm screw :40 mm
H15 SA6045	Ø 6 mm screw :45 mm
H15 S015	Peg Ø 6 mm Lg 15 mm
H29 005	5 Cup plugs

## PEAKED ANTI-LUXATION POLYETHYLENE INLAYS

Order ref.	Description for 28 & 32 PE inlays
H29 2246	Inlay for 22 mm head for cup size 46
H29 2848	Inlay for 28 mm head for cup size 48
H29 2850	Inlay for 28 mm head for cup size 50
H29 2852	Inlay for 28 mm head for cup size 52 & 54
H29 3252	Inlay for 32 mm head for cup size 52 & 54
H29 2856	Inlay for 28 mm head for cup size 56 & 58
H29 3256	Inlay for 32 mm head for cup size 56 & 58
H29 2860	Inlay for 28 mm head for cup size 60 & 63
H29 3260	Inlay for 32 mm head for cup size 60 & 63
H29 2866	Inlay for 28 mm head for cup size 66 & 69
H29 3266	Inlay for 32 mm head for cup size 66 & 69


## ALUMINA CERAMIC INLAY

Order ref.	Product description
H29 A2848	Alumina ceramic inlay for 28 mm head in cup size 48 mm
H29 A2850	No longer available. Only H29 A3250 is now used.
H29 A3250	Alumina ceramic inlay for 32 mm head in cup size 50
H29 A3252	Alumina ceramic inlay for 32 mm head in cup size 52 mm & 54
H29 A3256	Alumina ceramic inlay for 32 mm head in cup size 56 & 58 mm
H29 A3260	Alumina ceramic inlay for 32 mm head in cup size 60 & 63 mm
H29 A3266	Alumina ceramic inlay for 32 mm head in cup size 66 & 69 mm

## FLAT CONICAL PE INLAYS

H29 P2244	<i>Spherical</i> inlay for 22.2 mm head for cup size 44
H29 P2246	Inlay for 22.2 mm head for cup size 46
H29 P2848	Inlay for 28 mm head for cup size 48
H29 P2850	Inlay for 28 mm head for cup size 50
H29 P2852	Inlay for 28 mm head for cup size 52 & 54
H29 P2856	Inlay for 28 mm head for cup size 56 & 58
H29 P2860	Inlay for 28 mm head for cup size 60 & 63
H29 P2866	Inlay for 28 mm head for cup size 66 & 69

Rationalized and simplified instrument set to make facilitate operating and to ease with cleansing.



### IMPORTANT :

Prior to inserting the Alumina ceramic inlay check that the inside of the cup is clean. The use of the suction cup instrument is essential in order to ensure a correct fit.

To be able to remove the device without compromising the fixation of the inlay pull trigger to release suction vacuum.



**CE**  
**0499**

A SCANOS product  
☎ : +33.5.61.73.67.16 ☎ : +33.5.61.73.52.35  
e-mail : [biomed.biomaterials@wanadoo.fr](mailto:biomed.biomaterials@wanadoo.fr)  
visit us at : [www.scanos.com](http://www.scanos.com)

Manufactured by Evolutis S.A.S  
Avenue de la Libération - F- 42720 BRIENNON

REF H29PF 501 - 11/09