

SCANOS FEMORAL COMPONENT ORDER REFERENCES

HACTIV non cemented dysplastic HAC stem

Order ref.	Product description
H49 007	Hactiv non cemented dysplasia stem size 7
H49 008	Hactiv non cemented dysplasia stem size 8

HACTIV non cemented HAC stem

H49 009	Hactiv non collared HAC stem size 9
H49 010	Hactiv non collared HAC stem size 10
H49 011	Hactiv non collared HAC stem size 11
H49 012	Hactiv non collared HAC stem size 12
H49 013	Hactiv non collared HAC stem size 13
H49 014	Hactiv non collared HAC stem size 14
H49 015	Hactiv non collared HAC stem size 15
H49 016	Hactiv non collared HAC stem size 16
H49 018	Hactiv non collared HAC stem size 18
H49 020	Hactiv non collared HAC stem size 20

HACTIV non cemented HAC collared stem

H49 C009	Hactiv non cemented collared HAC stem size 9
H49 C010	Hactiv non cemented collared HAC stem size 10
H49 C011	Hactiv non cemented collared HAC stem size 11
H49 C012	Hactiv non cemented collared HAC stem size 12
H49 C013	Hactiv non cemented collared HAC stem size 13
H49 C014	Hactiv non cemented collared HAC stem size 14
H49 C015	Hactiv non cemented collared HAC stem size 15
H49 C016	Hactiv non cemented collared HAC stem size 16
H49 C018	Hactiv non cemented collared HAC stem size 18
H49 C020	Hactiv non cemented collared HAC stem size 20

HACTIV High Off-Set coxa vara stem

H49 L009	Hactiv high offset coxa vara HAC stem size 9
H49 L010	Hactiv high offset coxa vara HAC stem size 10
H49 L011	Hactiv high offset coxa vara HAC stem size 11
H49 L012	Hactiv high offset coxa vara HAC stem size 12
H49 L013	Hactiv high offset coxa vara HAC stem size 13
H49 L014	Hactiv high offset coxa vara HAC stem size 14
H49 L015	Hactiv high offset coxa vara HAC stem size 15
H49 L016	Hactiv high offset coxa vara HAC stem size 16

HACTIV High Off-Set collared coxa vara

H49 LC009	Hactiv high offset coxa vara collared HAC size 9
H49 LC010	Hactiv high offset coxa vara collared HAC size 10
H49 LC011	Hactiv high offset coxa vara collared HAC size 11
H49 LC012	Hactiv high offset coxa vara collared HAC size 12
H49 LC013	Hactiv high offset coxa vara collared HAC size 13
H49 LC014	Hactiv high offset coxa vara collared HAC size 14
H49 LC015	Hactiv high offset coxa vara collared HAC size 15
H49 LC016	Hactiv high offset coxa vara collared HAC size 16

HACTIV sss HAC stem (under development)

Order ref.	Product description
H49 005	Hactiv non cemented small stature stem size 5
H49 006	Hactiv non cemented small stature stem size 6

REACTIV coronal slotted revision HA coated stem

Order ref.	Product description	Length
H49 R012	REACTIV slotted stem without collar size 12	180 mm
H49 R014	REACTIV slotted stem without collar size 14	200 mm
H49 R016	REACTIV slotted stem without collar size 16	210 mm
H49 R018	REACTIV slotted stem without collar size 18	220 mm
H49 R020	REACTIV slotted stem without collar size 20	240 mm

REACTIV coronal slotted collared revision stem

Order ref.	Product description	Length
H49 RC012	REACTIV slotted stem with collar size 12	180 mm
H49 RC014	REACTIV slotted stem with collar size 14	200 mm
H49 RC016	REACTIV slotted stem with collar size 16	210 mm
H49 RC018	REACTIV slotted stem with collar size 18	220 mm
H49 RC020	REACTIV slotted stem with collar size 20	240 mm

REACTIV distal locking collared revision stem

Order ref.	Product description	Length
H49 RC112	REACTIV distal locking collared size 12	180 mm
H49 RC114	REACTIV distal locking collared size 14	200 mm
H49 RC116	REACTIV distal locking collared size 16	210 mm
H49 RC118	REACTIV distal locking collared size 18	220 mm
H49 RC120	REACTIV distal locking collared size 20	240 mm

CEMTIV stainless steel M30NW cemented stem

H49 S09	Cemtiv stainless steel cemented revision stem size 9
H49 S10	Cemtiv stainless steel cemented revision stem size 10
H49 S11	Cemtiv stainless steel cemented revision stem size 11
H49 S12	Cemtiv stainless steel cemented revision stem size 12
H49 S13	Cemtiv stainless steel cemented revision stem size 13
H49 S14	Cemtiv stainless steel cemented revision stem size 14
H49 S15	Cemtiv stainless steel cemented revision stem size 15
H49 S16	Cemtiv stainless steel cemented revision stem size 16

CEMTIV stainless steel collared revision stem

Order ref.	Product description	Length
H49 SRC012	Cemtiv revision stem with collar size 12	180 mm
H49 SRC014	Cemtiv revision stem with collar size 14	200 mm
H49 SRC016	Cemtiv revision stem with collar size 16	210 mm
H49 SRC018	Cemtiv revision stem with collar size 18	220 mm

NB: The REACTIV distal locking revision stem is the property of EVOLUTIS S.A.S and completes the product range. In the same manner the High Offset version of the HACTIV femoral component was expressly designed for X-Nov (France) and the anatomical Keratos version designed for OTHESIO (France).

All 3 can be supplied to you by Biomed-Scanos. All the other versions are the exclusive property of Biomed (France)